




EMSP

Environmental Management Science Program



Project #42370. Radioanalytical Chemistry for Automated Nuclear Waste Process Monitoring

Oleg Egorov, PNNL

Jay Grate, PNNL

Timothy DeVol, Clemson University



EMSP

Environmental Management Science Program



Radioanalytical Chemistry for Automated Nuclear Waste Process Monitoring

- Research program directed toward rapid, sensitive and selective determination of β - and α - emitters in LAW processing streams:
 - specific targets ^{90}Sr , ^{99}Tc and TRU
- Automated radiochemical measurements via integration of A) *sample preparation/treatment*, B) *rapid selective separation*, and C) *on-line radiometric detection* steps within a single functional analytical instrument
- Knowledge gaps exist relative to designing chemistries for such instruments so that analytes can be quantitatively and rapidly separated and analyzed in challenging matrixes



EMSP

Environmental Management Science Program



Project Objectives

- Introduce automated radiochemistry to the field of process monitoring of LAW matrices
- Develop fundamental understanding of the chemical processes that must occur within the analyzer instrument (i.e., matrix modification/speciation control and selective separation chemistries)
- Investigate matrix modification chemistries for chemically complex, variable matrices
- Investigate separation chemistries suitable for rapid selective separation of target analytes
- Investigate a new flow-through radiometric detector concept based on solid state diode detectors for quantification of TRU in process solutions
- Investigate new data acquisition and signal processing techniques to enhance performance of scintillation and diode detectors for use in process monitoring
- Characterize candidate chemistries and instrument configurations with regard to detection limits, sensitivity, selectivity, analyses time
- Define the limitations and the potential of the proposed monitoring approach. New knowledge and expertise will provide information for designing effective instruments
- Educate students in modern radiochemistry and DOE cleanup challenges



EMSP

Environmental Management Science Program




Process Monitoring Challenges and Reality

- Generic Challenges:
 - Complex, high-ionic strength, caustic brine sample matrixes
 - Undefined or varying chemical speciation
 - Low mass concentrations, need for isotopic information
 - Radioactive and stable matrix interferences
 - Detection limits, robustness, sample turnaround time
- Challenges/requirements limit the choice of analytical approaches
- **Reality:** “Radiochemistry has always been and still is a crucial tool in the field of radionuclide determination, ...particularly in the case of alpha and beta emitters”.
- Current baseline: costly and lengthy laboratory analyses





EMSP

Environmental Management Science Program



Automated Radiochemistry: Approach




■ Automated radiochemical measurements via integration of the following steps within a single functional instrument:

- sample preparation/treatment,
- rapid selective separation,
- on-line radiometric detection



EMSP

Environmental Management Science Program




Automated Radiochemistry: Research Focus

- Sample modification/ speciation control chemistries for ^{99}Tc , ^{90}Sr (TRU?) (acidification, digestion, oxidation, speciation control)
- Separation chemistries for ^{99}Tc and ^{90}Sr (hopefully not for TRU)
 - sorbent material selection and characterization (uptake selectivity, kinetics, elution chemistries, stability etc.
 - development of rapid, automated separation schemes
- On-line radiometric detection
 - solid state diode detection for in-situ alpha detection/spectroscopy
 - light detection, detector design and signal processing for beta scintillation detection (Clemson)
- Analytical and radiochemistry of the LAW waste matrices



EMSP

Environmental Management Science Program




Relevance

- ^{90}Sr , ^{99}Tc and TRU process measurements will be required at WTP at Hanford and Savannah River sites
 - Better process control, plant throughput and regulatory compliance
- Lack of analytical methodologies for these radionuclides
- Baseline is not significantly different from PUREX plant days
- Plants are being constructed: this research is timely
- This analytical methodology has dual use:
 - on-line/at-line measurements
 - automated tool for the analytical laboratory
- Significant cost savings over operational cycle of the waste treatment plants



EMSP


Environmental Management Science Program



Automated Radiochemistry: Challenges

Tank AN-102 -Radiochemical interference issue in the automated Tc analysis.


Inventory in	uCi/mL
Gross alpha:	0.18
Gross beta:	590
Tc-99	8.23(ppm)
(at least 60% as non-pertechnetate)	
Cs-137	335
Cs-135	3.02
Sr-90	84
Co-60	0.12
Pu-239/240	0.01
Am-241	0.15
Unknown Activity	158






EMSP

Environmental Management Science Program



Automated Radiochemistry Example: Tc Process Monitoring




- **Flow-through scintillation detection of the beta-emitting ⁹⁹Tc**
(radiometric detection is sensitive to meet DL with short counting times)
 - Detection using glass scintillator flow cell
- **Automated, rapid radiochemical separation of the ⁹⁹Tc(VII) from sample matrix and radioactive interferences (e.g. ⁹⁰Sr/⁹⁰Y, ¹³⁷Cs, etc.)**
(separation is required to measure ⁹⁹Tc in the presence of other radionuclides)
 - Small column separation using anion exchange sorbent
- **Automated, rapid sample treatment/oxidation to control Tc speciation**
(oxidation is required to convert all ⁹⁹Tc to ⁹⁹Tc(VII) for total ⁹⁹Tc analysis)
 - Sample acidification and microwave assisted oxidation using persulfate
- **Automated, digital fluid handling for sample/solution delivery and system integration**
 - Computer-controlled valves and syringe pumps





EMSP

Environmental Management Science Program




Prototype Instrument for Hanford WTP






EMSP

Environmental Management Science Program



Microwave Assisted Sample Oxidation



- Unknown speciation
- Reported difficulties for converting Tc to Tc(VII)
- Need for rapid and reliable methodology
- Must be practical for use in plant settings
- MW –assisted format using acidification and oxidation



EMSP

Environmental Management Science Program



Microwave Assisted Sample Oxidation for Tc Speciation Control

1) Sample acidification using nitric acid (nitrite removal)

- $\text{Al}(\text{OH})_3$ precipitation/dissolution;
- Digested matrix must be compatible with subsequent separation approach

2) Tc oxidation using $\text{Na}_2\text{S}_2\text{O}_8$ as the oxidizing reagent

- Complete destruction of organic complexants is not necessary
- **Ag catalyst is typically used but not necessary**
- No Tc losses during digestion
- Rapid, reliable, automated process


Advantages of peroxodisulfate:

- strong oxidant ($E^0=2.01$ V)
- good shelf life
- excess reagent is readily destroyed by boiling





EMSP

Environmental Management Science Program



Microwave-Assisted Oxidation of Non-Per technetate

Envelope C AN-107 waste: ~ 60% of non-per technetate
Chemical speciation is unknown




• Rapid, reliable Tc oxidation is possible using MW-assisted $\text{Na}_2\text{S}_2\text{O}_8$ treatment




EMSP

Environmental Management Science Program




Identification of Initial Interferent in Automated Radiochemical Tc Measurement

- liquid scintillation spectrometry of the separated Tc fractions indicated good radiochemical purity for the analysis of radiochemically complex samples



Additional verification:


LEPS
Gamma
ICP MS






EMSP

Environmental Management Science Program



Automated Tc Measurements in Hanford Waste




- High non-pertechetate content
- Chemically & radiochemically most challenging matrices
- Also verified for Envelope A, B, and C wastes

- Samples were batch contacted with SL639 Resin
- Better than 10% agreement for both matrixes
- No bias




EMSP

Environmental Management Science Program



Tc Process Monitor Instrument Testing for WTP



Sample	Accuracy (%)	Precision (%RSD)	N
60XAZ102	4.4	6.2	45
50XAZ102	2.0	7.6	50
40XAZ102	1.2	6.6	34
30XAZ102	2.0	3.8	31

•Analytical approach/instrument selected for Hanford WTP in peer reviewed, competitive selection process!!!




EMSP

Environmental Management Science Program



Direct Detection of Alpha Emitters In Solution Using PIPS Diodes


- Low noise and good resolution
- Operation in ambient light
- Passivated surface w/o metal contacts
 - Rugged
 - Chemically resistant
 - Surface can be touched, cleaned--*modified*






EMSP

Environmental Management Science Program




Direct Analysis of Alpha Emitters In Solution

Activity vs. Distance*



* Distance between solution surface and detector face. 125 nCi/ml ^{241}Am activity.

Alpha Energy "Spectra"†




† 125 nCi/ml, 10 hr counts, 0.5 cm standoff from surface.




EMSP

Environmental Management Science Program



Direct AEA in Liquids by Spectral Deconvolution



Isotope	Solution Activity ^a	Calculated Activity ^b
²⁴⁴ Cm	25.0	24.8
²⁴¹ Am	18.8	20.2
²³⁹ Pu	62.5	60.9
²³³ U	18.8	19.3

- a) nCi/mL, mixture of standards
- b) nCi/mL, spectral deconvolution, least squares optimization

20% ²⁴⁴Cm, 15% ²⁴¹Am, 50% ²³⁹Pu, 15% ²³³U
Total solution activity 125 nCi/mL



EMSP

Environmental Management Science Program




Analysis of Actual Samples

Sample Matrix	Density (g/mL)	PIPS Assay ^a (Bq/ml)	Lab Assay ^b (Bq/ml)	Difference (%)
Water ^c	1	4.6	4.0	-13.9
Water ^c	1	392.7	403.0	2.6
Hanford Tank AZ-101 ^{c,d}	1.22	17.1	15.7	-9.2
Hanford Tank AN-107 ^d	1.249	70.4	82.1	14.2
Hanford Tank AZ-102 ^d	1.260	120.0	154.9	22.5
Savanna River Tank 23 ^{c,d}	1.031	817.6	745.0	-9.7
Savanna River Tank 33 ^d	1.171	6.6	8.6	22.7
Savanna River Tank 43 ^d	1.204	224.8	205.5	-9.4

- a) Based upon ²⁴¹Am, 10 min count, efficiency corrected
- b) Precipitation plating and AES
- c) Spiked with ²⁴¹Am
- d) supernate

- Rapid total alpha assay
- Limited sensitivity ~5 Bq/mL
- Spectral information....





EMSP

Environmental Management Science Program




Implementation

- Hanford WTP:
 - On-line Tc monitoring is required for the operation of the Tc removal columns
 - Our technology was selected for the Hanford WTP
 - Sr and TRU monitoring: automated laboratory measurements on grab samples (so far no interest at Hanford)
- Salt processing facility at SRS
 - Need for rapid Sr and TRU analysis in the lab to support process control
 - Indication that such technology may be a preferred method for this need
 - Goal is to pursue prototype instrumentation development with SRS using **new scientific knowledge** developed under this EMSP program



EMSP

Environmental Management Science Program



Current Focus & Conclusions

- Alternative separation chemistries for Sr and Tc
 - Simplify separation protocols, Ba/Sr separation issues
- Diode detection and signal processing
- Radiochemistry represent a viable approach for process monitoring if successfully automated
- Tc monitor instrumentation/ approach was selected for use at the Hanford WTP
- Work with WTP contractors would have not been possible without the science base and knowledge developed under an on going EMSP program



EMSP

Environmental Management Science Program



Acknowledgements

PNNL:

Matt O'Hara, Shane Addleman, Bruce McNamara,
Sergei Sinkov

CHG: Michael Johnson

BNI: Reid Peterson, David Sheerwood